

niccons

ACID TEST

Cod: TESTACID



DESCRIPTION

The acid test works with all lubricants and quickly and economically allows to detect the acidity of the lubricant inside the system. It is easy to use and quickly detects acid without having to wait weeks for a laboratory analysis, so it allow to save time and money. It's especially sensitive to small acid particles. Package valid for 2 tests. It includes: 3 measuring pipettes 5 ml, 1 bottle with 10 ml of determiner and 2 bottles with 1 ml of activator.

INSTRUCTIONS FOR USE

- Fill one of the two containers (containing the activator) with 10 ml of the oil to be examined.
- Using the pipette, slowly add and shake 0,5 ml of determiner at a time.
- Stop when the solution became purple.

The acidity will result:

- Low if we add only 0,5 ml of determiner (purple coloration reached after one step);
- Average if we add 1 ml of determiner (purple coloration reached after two steps);
- High if we add 1,5 ml to 5 ml of determiner (purple coloration reached from 3 to 10 steps). If the purple coloration is reached after putting in the activator container more than 5 ml of determiner ($> 0,5 \times 10$ steps) we need to proceed with the substitution of the oil and the complete washing of the system.

SAFETY INFORMATION

May be fatal or cause blindness if swallowed. May cause eye and skin irritation. Keep away from heat, sparks and flames. Use with adequate ventilation. Wash after handling. Store in a cool, well ventilated area, away from all food. Keep container upright and tightly closed. Keep out of reach of children. Do not reuse empty container. Professional use only.

NICCONS Italy srl

Via dell'Economia, 85 - 36100 Vicenza (VI)

Tel.: +39 0444 289383 - Email: info@niccons.com - Web: www.niccons.com